



Revolutionary Formulation of the state of th **ROOT CANAL BIOCERAMIC SEALER**

Breakthrough Technology with Strontium Silicate

Global Patent: PCT/CN2019/097075

FDA Clearance: K212983

Registration/FEI Number: 3021971019 Owner/Operator Number: 10084370

US PATENT NO.17/266.102 EU PATENT NO.19847130.2



Scan QR Codes & Get free sample















huangmeijuan@crootmed.com







Product Description

- C-Root SP Injectable Root Canal Bioceramic Sealer is a premixed, injectable silicate bio-ceramic paste material which solidifies in the presence of water. It is used to sealing the various complexities within the root canal system during obturation and promote mineralization. C-Root SP is a mixture containing strontium silicate, calcium, silicon, strontium and magnesium ions, which harden in the presence of water and promotes healing and dentin mineralization; It is bactericidal due to its high alkalinity prior to it final setting.
- C-Root SP is a newly patent awarded Bioceramic Root Canal Sealer. It is injectable white hard cement, which is designed for safe and effective permanent root canal filling and sealing.

Composition: Strontium Silicate, Zirconia, Calcium Phosphate and Calcium Hydroxide

Advantages: Excellent biocompatibility. Easy to use. Lateral canal filling ability. Excellent radiopacity. Mild postoperative reaction. Rapid periapical bone lesion healing. No shrinkage during setting.

Used for: Root Canal Sealing and Permanent Filling.

Indications:Permanent filling of root canal after it has been cleaned and disinfected properly.

It is suitable for lateral condensation/compaction technique, single cone obturation technique and vertical condensation techniques



Device Listing:D466741



Strontium Silicate Bioceramics Experience the extraordinary

C-ROOT SP INJECTABLE ROOT CANAL BIOCERAMIC SEALER Biocompatibility and High X-ray radio opacity 3D sealing and filling Osteogenesis Hydrophilicity Easy operation Chemically bonds to dentin React to form hydroxyapatite High antibacterial (pH > 12) No shrinkage Remineralization in dental **Affordability** Ideal working time tubules Lateral canal filled Rapid periapical healing Single cone technique Average fill 15-25 root Mild Post-Op Reaction **Excellent flow characteristics** canals per 2g





Product Properties

 $\cap 1$

GOOD SEALING, REMINERALIZATION IN DENTAL TUBULES TO GENERATE STRONTIUM HYDROXYAPATITE

- ➤ C-Root SP particles are relatively fine, which can be injected into the root canal, flow into the various anatomical complexities (including small lateral root canal, fins, cul-de-sacs). It can also penetrate the dental tubules (2-5 microns in diameter) for remineralization to form hydroxyapatite with strontium substituted, which effectively seals the canal.
- C-Root SP bioceramic premixed paste, is easy to operate and can be filled by single cone technique, which can save time and reduce the risk of root fracture. The price is reasonable and acceptable, it has great application value.





03

ANTIBACTERIAL

Antibacterial efficacy of 100% of enterococcus feacalis and 99% of staphylococcus Aureus within 2 minutes of exposure. C-Root SP contains calcium hydroxide which is a potent bactericidal agent commonly used in root canal medicaments which enhances its bactericidal properties.

时间节点	杀灭对象	杭菌效果					
		×1		×10		×100	
2min	试验样品组重复 1						
	金黄色葡萄球菌	99%	100%	99%	100%	100%	100%
	黄肠球菌	100%	100%	100%	100%	100%	100%
	试验样品组重复 2						
	金黄色葡萄球菌	100%	100%	100%	100%	100%	100%
	黄肠球菌	100%	100%	100%	100%	100%	100%
	试验样品组重复 3						
	金黄色葡萄球菌	100%	100%	100%	100%	100%	100%
	黄肠球菌	100%	100%	100%	100%	100%	100%
	试验样品组重复 4						
	金黄色葡萄球菌	100%	100%	100%	100%	100%	100%
	类肠球菌	100%	100%	100%	100%	100%	100%


02

GOOD BIOCOMPATIBILITY

C-Root SP has been tested for biocompatibility by The Dental medical device testing center, School of Stomatology, Peking University for up to 2 years and reported to be non-toxic, non irritating, non carcinogenic and non mutagenic with osteogenenic potential and mild post operative reactions.


Endodontic Usage Test (REP/PDS64-2019-0014-GG-BGE) showed:

At the time period of 90 days, among the 10 test specimens (C-Root SP),10 specimens (C-Root SP) showed no inflammation.



Among the 5 control specimens (other sealer), 4 specimens showed no inflammation, 1 specimens of periapical tissue showed moderate inflammation (See the figure below).






90d- Control- 1 (25X)

90d- Control- 1 (100X)

04

INDUCE OSTEOGENESIS AND ROOT END FORMATION

C-Root SP is not easy to extrude compared to other sealers! Small amount of extrusion can induce osteogenesis and root end formation. Good apical sealing and osteogenic potential from the strontium silicate which promotes healing for long term success.



CLEARED BY FDA IN K212983

Pre-Clinical Research Report

In vitro cytocompatibility and osteogenic potential of calcium silicate-based dental cements in a root canal-filling model

BELING C-ROOT DENTAL MEDICAL DEVICES CO

Yan Jing^{1,2}, Ting Gong², Chenggang Duan¹, Huan Wang¹, Chengfei Zhang² and Prasanna Neelakantan²

Abstract

Objective: To investigate the in vitro cytocompatibility and osteogenic potential of an experimental calcium silicate-based cement and the inflammatory response in human periodontal ament stem cells (hPDLSCs).

Methods: Cellular responses, osteogenic-related gene expression, and the production of inflammatory cytokines including interleukin (III). As and II-8 were studied in hPDLSCs exposed to the experimental root canal-filling material C-Root, the commercial tricalcium silicate-based material BioRoot RCS, and the epoxy resin-based material AH Plus. Differences were analyzed using one-way ANOVA with Bonferroni-adjusted pairwise comparison.



WHAT IS THE IDEAL ROOT CANAL BIOCERAMIC SEALER

Flow Setting time Sealing 3D filling X-ray radio opacity No staining Bactericidal Biocompatibility No resorbable Mild Post-OP Osteogenesis Remineralization in dental tubules Reaction







Located in Mapo Science and Technology Park, Shunyi District, Beijing. Beijing C-Root Dental Medical Devices Co. LTD is a Zhongguancun high-tech enterprise with R & D, manufacture and sale of dental root canal treatment devices, which has achieved breakthrough development.

C-Root Dental Company has built a GMP quality management system and a purification workshop of 600 square meters; Two types of medical devices, EDTA root canal lubricating gel and Calcium hydroxide root canal medication paste, have been highly recommended by general dentists and Endodontic specialists.

C-Root SP Injectable Root Canal Bioceramic Sealer, has been awarded the registration certificate of CFDA (gxzz 20213170769) of China and FDA certificate (K212983) of the United States of America . C-Root SP has global

patent protection for strontium silicate paste. The excellent properties of C-Root SP ensure high and long term success, and the affordable price benefits dentists and patients.





C-Root NiTi Files

Flexible, Fast Root Canal Preparation

C-Root NiTi Files .04:

015/08 19mm 2.5Necm 500RPM 015/03 21mm 1.5Necm 300RPM 020/04 21mm 2.5Necm 500RPM 025/04 21mm 2.5Necm 500RPM

C-Root NiTi Files .06:

015/08 19mm 2.5Necm 500RPM 015/03 21mm 1.5Necm 300RPM 020/05 21mm 2.5Necm 500RPM 025/06 21mm 2.5Necm 500RPM

Device Listing: D466740





EDTA ROOT CANALLUBRICATING GEL

Lubrication and file fracture prevention.



Composition: 10% urea peroxide + 17% EDTA + glycerol + polyethylene glycol + carboxymethyl cellulose CMC + sodium hyaluronate


DISINFECTION PASTE

The recognized most effective root canal



- √ Good flow
 √ X-ray radio opacity
- No setting V Easy to take out
- Easy to put into root tip 1/3 part
- Complete sterilization, less pain during diagnosis and less postoperative reaction


WHAT KIND OF FILLING DO YOU WANT?





Chemical Remineralization C-Root SP of Dentinal Tubules Bongding











By Dr. Sun Yi , C-Root SP





Is It Still Your Choice?















& Get free sample

Contact: crootsp (+86) 18601290788

www.crootmed.com

Beijing C-Root Dental Medical Devices Co., LTD.

Registrant / manufacturer: Beijing C-Root Dental Medical Devices Co., LTD Add: Room 301, Building 10, No 12, Juyuan Middle Road, Mapo Town, Shunyi , Beijing, China 101300 Tel: (+86) 13371625696 / (+86) 18601290788

Email: wubingmin@crootmed.com / huangmeijuan@crootmed.com